
**ADVICES ON CHECKS AND CLEANING
REMOVABLE DENTAL PROSTHESIS
VIVA FLEX PARTIALS**



AUTHORIZED SERVICE

Technical assistance in the event there is damage to the dental prosthesis can be carried out by the professional who applied it or by the staff of the manufacturer Dental Flex Italia S.n.c..

Checks

Periodically it is necessary to control the prosthesis. To do this you must:

- Wash the prosthesis at least twice a day, especially after meals.
- Brushing the prosthesis VIVA FLEX PARTIALS with a normal toothbrush under running water.
- It is advisable to carry out periodic checks on the device at least every six (6) months, unless introducing subtle acts of God.
- If the prosthesis monsters plaque buildup or smudges, due to a NOT correct oral hygiene. It is important for the doctor to send back the prosthesis laboratory that manufactured the same; they will be able to perform cleaning or mechanical.

Routine cleaning

At least two (2) times a day, conduct the normal oral hygiene by washing the teeth. Use conventional toothbrushes and toothpaste.

Take care of your new prosthesis VIVA FLEX PARTIALS!

Information on the manufacturer

Name Legal Representative	Gianfranco
Surname Legal Representative	Guerra
Enterprise	Dental Flex Italia S.n.c.
Company headquarters: Street and City	Via Montevergine, 62 – Ospedaletto d’Alpinolo (AV) – 83014 - Italia
VAT number	02570240644
Product name	Removable dental prosthesis VIVA FLEX PARTIALS
Tel. Office	+39 0825 - 691237
Fax.	+39 0825 - 691237
Email	gianfranco@dentalflexitalia.eu
Website	www.dentalflexitalia.eu
Product link (if available)	